

希釈冷凍機 ポート 1 温度計チェックシート No. 1

-ishisaki/tx/dilution/tejun/port1_check1.pdf

日付：

記入者

コメント

Function	Pin No	抵抗[Ω]	
IGH Sorb sensor I (455 Ω)	1-18		
IGH Sorb sensor V (455 Ω)	2-17		
IGH 1K pot sensor I (2378 Ω)	3-18		
IGH 1K pot sensor V (2388 Ω)	4-17		
IGH Mixing chamber I (2434 Ω)	9-18		
IGH Mixing chamber V (2434 Ω)	10-17		
IGH - GND	18		
AVS Still I (2418 Ω)	5-16		
AVS Still V (2417 Ω)	6-15		
AVS Cold plate I (2440 Ω)	7-16		
AVS Cold plate V (2439 Ω)	8-15		
AVS Cold plate + (178 Ω)	7-8		
AVS common - (177 Ω)	15-16		
AVS - GND	16		
AVS Mixing chamber I (2474 Ω)	11-14		
AVS Mixing chamber V (2476 Ω)	12-13		
AVS Mixing chamber + (260 Ω)	11-12		
AVS Mixing chamber - (261 Ω)	13-14		
AVS Mixing chamber - GND	11		
Sorb heater (72 Ω)	19-20		
Sorb heater - GND	20		
Still heater (539 Ω)	21-22		
Still heater - GND	22		
M/C heater (555 Ω)	23-24		
M/C heater - GND	24		